

A5047

Leader in Biomolecular Solutions for Life Science



# HSDL2 Rabbit pAb

Catalog No.: A5047

## Basic Information

---

### Observed MW

46kDa

### Calculated MW

45kDa

### Category

Mouse Monoclonal Antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

---

Predicted to enable oxidoreductase activity. Located in mitochondrion and peroxisome.

## Recommended Dilutions

---

WB 1:500 - 1:1000

## Immunogen Information

---

### Gene ID

84263

### Swiss Prot

Q6YN16

### Immunogen

A synthetic Peptide of human HSDL2

### Synonyms

C9orf99; SDR13C1; HSDL2

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide, pH7.3.

## Validation Data

---



Western blot analysis of extracts of 293T cells, using HSDL2 antibody (A5047).  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.