

ACCS Rabbit pAb

Catalog No.: A5062

Basic Information

Observed MW

54kDa

Calculated MW

57kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables identical protein binding activity. Predicted to be involved in biosynthetic process.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

84680

Swiss Prot

Q96QU6

Immunogen

A synthetic Peptide of human ACCS

Synonyms

ACS; PHACS; ACCS

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of MDA-MB-453 cells, using ACCS antibody (A5062).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.