

# HSF2 Rabbit mAb

Catalog No.: A5090 **Recombinant**

## Basic Information

---

### Observed MW

### Calculated MW

60kDa

### Category

Primary antibody

### Applications

0

### Cross-Reactivity

## Background

---

The protein encoded by this gene belongs to the HSF family of transcription factors that bind specifically to the heat-shock promoter element and activate transcription. Heat shock transcription factors activate heat-shock response genes under conditions of heat or other stresses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

## Recommended Dilutions

---

0

## Immunogen Information

---

### Gene ID

3298

### Swiss Prot

Q03933

### Immunogen

A synthesized peptide derived from human HSF2

### Synonyms

HSF 2; HSTF 2; HSF2

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.