

A5090

Leader in Biomolecular Solutions for Life Science



HSF2 Rabbit mAb

Catalog No.: A5090

Recombinant

Basic Information

Observed MW

Calculated MW

60kDa

Category

Mouse Monoclonal Antibody

Applications

Cross-Reactivity

Background

The protein encoded by this gene belongs to the HSF family of transcription factors that bind specifically to the heat-shock promoter element and activate transcription. Heat shock transcription factors activate heat-shock response genes under conditions of heat or other stresses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

Immunogen Information

Gene ID

3298

Swiss Prot

Q03933

Immunogen

A synthesized peptide derived from human HSF2

Synonyms

HSF 2; HSTF 2; HSF2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.