

ACAP3 Rabbit pAb

Catalog No.: A5123

Basic Information

Observed MW

Refer to figures

Calculated MW

92kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Predicted to enable GTPase activator activity and metal ion binding activity. Predicted to be involved in regulation of catalytic activity. Predicted to act upstream of or within neuron migration and regulation of neuron projection development. Predicted to be located in growth cone.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

116983

Swiss Prot

Q96P50

Immunogen

A synthetic Peptide of human ACAP3

Synonyms

CENTB5; ACAP3

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.