

A5320

Leader in Biomolecular Solutions for Life Science



NOV Rabbit pAb

Catalog No.: A5320

Basic Information

Observed MW

Calculated MW

39kDa

Category

Polyclonal Antibody

Applications

IHC-P,ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene is a small secreted cysteine-rich protein and a member of the CCN family of regulatory proteins. CNN family proteins associate with the extracellular matrix and play an important role in cardiovascular and skeletal development, fibrosis and cancer development.

Recommended Dilutions

IHC-P 1:50 - 1:100

Immunogen Information

Gene ID

4856

Swiss Prot

P48745

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 32-357 of human NOV (NP_002505.1).

Synonyms

NOV; NOVh; IBP-9; IGFBP9; IGFBP-9

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Immunohistochemistry analysis of paraffin-embedded Human prostate using NOV Rabbit pAb (A5320) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.