

A5575

Leader in Biomolecular Solutions for Life Science



RABL2A Guinea pig pAb

Catalog No.: A5575

Basic Information

Observed MW

Refer to figures

Calculated MW

26kDa/27kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Enables GTP binding activity. Involved in photoreceptor cell morphogenesis. Acts upstream of or within cilium assembly; flagellated sperm motility; and single fertilization. Located in sperm midpiece. Part of intraciliary transport particle B. Is expressed in several structures, including eye; genitourinary system; gut; nervous system; and respiratory system. Orthologous to several human genes including RABL2B (RAB, member of RAS oncogene family like 2B).

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

68708

Swiss Prot

Immunogen

Recombinant Protein of mouse RABL2A

Synonyms

Rabl2a; 1110031N17Rik

Contact



www.abclonal.com

Product Information

Source

Guinea pig

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.