

GRWD1 Rabbit pAb

Catalog No.: A6002

Basic Information

Observed MW

49kDa

Calculated MW

49kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a glutamate-rich protein that contains five WD-repeat motifs. The encoded protein may play a critical role in ribosome biogenesis and may also play a role in histone methylation through interactions with CUL4-DDB1 ubiquitin E3 ligase.

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

83743

Swiss Prot

Q9BQ67

Immunogen

A synthetic peptide of human GRWD1

Synonyms

CDW4; GRWD; RRB1; WDR28; GRWD1

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of 293 cells, using GRWD1 antibody (A6002).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.