

GNL2 Rabbit pAb

Catalog No.: A6805

Basic Information

Observed MW
90kDa

Calculated MW
84kDa

Category
Primary antibody

Applications
WB

Cross-Reactivity
Human

Background

Enables RNA binding activity. Predicted to be involved in ribosome biogenesis. Located in nucleolus.

Recommended Dilutions

WB 1:200 - 1:1000

Immunogen Information

Gene ID	Swiss Prot
29889	Q13823

Immunogen
A synthetic peptide of human GNL2

Synonyms
NGP1; Nog2; Nug2; Ngp-1; HUMAUANTIG; GNL2

Contact


 | www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage
Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of 293T cells, using GNL2 antibody (A6805).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.