

A7907

Leader in Biomolecular Solutions for Life Science



EG5/KIF11 Rabbit pAb

Catalog No.: A7907 **2 Publications**

Basic Information

Observed MW

130kDa

Calculated MW

119kDa

Category

Polyclonal Antibody

Applications

WB,IF/ICC,ELISA

Cross-Reactivity

Human

Background

This gene encodes a motor protein that belongs to the kinesin-like protein family. Members of this protein family are known to be involved in various kinds of spindle dynamics. The function of this gene product includes chromosome positioning, centrosome separation and establishing a bipolar spindle during cell mitosis.

Recommended Dilutions

WB	1:500 - 1:1000
IF/ICC	1:50 - 1:100

Immunogen Information

Gene ID

3832

Swiss Prot

P52732

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 797-1056 of human EG5/EG5/KIF11 (NP_004514.2).

Synonyms

EG5; HKSP; KNSL1; MCLMR; TRIP5; EG5/KIF11

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using EG5/KIF11 Rabbit pAb (A7907) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.



Immunofluorescence analysis of U2OS cells using EG5/KIF11 Rabbit pAb (A7907) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.