

A8636

Leader in Biomolecular Solutions for Life Science



# YY1AP1 Rabbit pAb

Catalog No.: A8636

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

88kDa

### Category

Polyclonal Antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

---

Predicted to enable transcription coregulator activity. Involved in cell differentiation; cell population proliferation; and regulation of cell cycle. Located in fibrillar center and nucleoplasm. Colocalizes with Ino80 complex.

## Recommended Dilutions

---

WB 1:500 - 1:1000

## Immunogen Information

---

### Gene ID

55249

### Swiss Prot

Q9H869

### Immunogen

A synthetic Peptide of human YY1AP1

### Synonyms

GRNG; HCCA1; HCCA2; YY1AP; YY1AP1

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide, pH7.3.