

# ACAP2 Rabbit pAb

Catalog No.: A8955

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

88kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human, Mouse

## Background

---

Enables GTPase activator activity. Acts upstream of or within actin filament-based process. Located in ruffle.

## Recommended Dilutions

---

WB 1:500 - 1:1000

## Immunogen Information

---

### Gene ID

23527

### Swiss Prot

Q15057

### Immunogen

A synthetic Peptide of human ACAP2

### Synonyms

CENTB2; CNT-B2; ACAP2

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.