

AFG3L2 Rabbit pAb

Catalog No.: A9020

Basic Information

Observed MW

Refer to figures

Calculated MW

89kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a protein localized in mitochondria and closely related to paraplegin. The paraplegin gene is responsible for an autosomal recessive form of hereditary spastic paraplegia. This gene is a candidate gene for other hereditary spastic paraplegias or neurodegenerative disorders.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

10939

Swiss Prot

Q9Y4W6

Immunogen

A synthetic Peptide of human AFG3L2

Synonyms

OPA12; SCA28; SPAX5; AFG3L2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.