

GIPC1 Rabbit pAb

Catalog No.: A9042

Basic Information

Observed MW

Refer to figures

Calculated MW

36kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

GIPC1 is a scaffolding protein that regulates cell surface receptor expression and trafficking (Lee et al., 2008 [PubMed 18775991]).

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

10755

Swiss Prot

O14908

Immunogen

A synthetic Peptide of human GIPC1

Synonyms

NIP; GIPC; IIP-1; OPDM2; TIP-2; SEMCAP; C19orf3; Hs.6454; GLUT1CBP; RGS19IP1; SYNECTIN; SYNECTIIN; GIPC1

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.