

A9223

Leader in Biomolecular Solutions for Life Science



ALKBH5 Rabbit pAb

Catalog No.: A9223

Basic Information

Observed MW

Refer to figures

Calculated MW

44kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables mRNA N6-methyladenosine dioxygenase activity. Involved in RNA metabolic process; mRNA export from nucleus; and response to hypoxia. Located in Golgi apparatus; cytosol; and nuclear speck.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

54890

Swiss Prot

Q6P6C2

Immunogen

A synthetic Peptide of human ALKBH5

Synonyms

ABH5; OFOXD; OFOXD1; ALKBH5

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.