

# CTDSPL2 Rabbit pAb

Catalog No.: A9275

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

53kDa

### Category

Primary antibody

### Applications

WB,FC

### Cross-Reactivity

Human, Mouse

## Background

---

Enables RNA polymerase II CTD heptapeptide repeat phosphatase activity. Predicted to act upstream of or within negative regulation of BMP signaling pathway; positive regulation of protein export from nucleus; and protein dephosphorylation. Located in nucleoplasm.

## Recommended Dilutions

---

WB	1:500 - 1:1000
FC	1:20 - 1:50

## Immunogen Information

---

<b>Gene ID</b>	<b>Swiss Prot</b>
51496	Q05D32

### Immunogen

A synthetic Peptide of human CTDSPL2

### Synonyms

SCP4; HSPC058; HSPC129; CTDSPL2

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

<b>Source</b>	<b>Isotype</b>	<b>Purification</b>
Rabbit	IgG	Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide,pH7.3.