

# DLL4 Rabbit pAb

Catalog No.: A9471

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

75kDa

### Category

Primary antibody

### Applications

WB,FC

### Cross-Reactivity

Human

## Background

---

This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain.

## Recommended Dilutions

---

<b>WB</b>	1:500 - 1:1000
<b>FC</b>	1:20 - 1:50

## Immunogen Information

---

<b>Gene ID</b>	<b>Swiss Prot</b>
54567	Q9NR61

### Immunogen

A synthetic Peptide of human DLL4

### Synonyms

AOS6; delta4; hdelta2; DLL4

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

<b>Source</b>	<b>Isotype</b>	<b>Purification</b>
Rabbit	IgG	Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide,pH7.3.