

CHST2 Rabbit pAb

Catalog No.: A9746

Basic Information

Observed MW

Refer to figures

Calculated MW

58kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

This locus encodes a sulfotransferase protein. The encoded enzyme catalyzes the sulfation of a nonreducing N-acetylglucosamine residue, and may play a role in biosynthesis of 6-sulfosialyl Lewis X antigen.

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

9435

Swiss Prot

Q9Y4C5

Immunogen

A synthetic Peptide of human CHST2

Synonyms

C6ST; GST2; GST-2; Gn6ST-1; HEL-S-75; glcNAc6ST-1; CHST2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.