

A9763

Leader in Biomolecular Solutions for Life Science



# CDK3 Rabbit pAb

Catalog No.: A9763

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

35kDa

### Category

Polyclonal Antibody

### Applications

WB,FC

### Cross-Reactivity

Human

## Background

---

This gene encodes a member of the cyclin-dependent protein kinase family. The protein promotes entry into S phase, in part by activating members of the E2F family of transcription factors. The protein also associates with cyclin C and phosphorylates the retinoblastoma 1 protein to promote exit from G0.

## Recommended Dilutions

---

WB	1:500 - 1:1000
FC	1:20 - 1:50

## Immunogen Information

---

<b>Gene ID</b>	<b>Swiss Prot</b>
1018	Q00526

### Immunogen

A synthetic peptide of human CDK3

### Synonyms

CDK3

## Contact

---

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

<b>Source</b>	<b>Isotype</b>	<b>Purification</b>
Rabbit	IgG	Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide,pH7.3.