

AP0042

Leader in Biomolecular Solutions for Life Science



Phospho-Insulin Receptor-Y1361 Rabbit pAb

Catalog No.: AP0042

Basic Information

Observed MW

110kDa

Calculated MW

156kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the receptor tyrosine kinase family of proteins. The encoded preproprotein is proteolytically processed to generate alpha and beta subunits that form a heterotetrameric receptor. Binding of insulin or other ligands to this receptor activates the insulin signaling pathway, which regulates glucose uptake and release, as well as the synthesis and storage of carbohydrates, lipids and protein. Mutations in this gene underlie the inherited severe insulin resistance syndromes including type A insulin resistance syndrome, Donohue syndrome and Rabson-Mendenhall syndrome. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3643

Swiss Prot

P06213

Immunogen

A synthetic phosphorylated peptide around Y1361 of human INSR (NP_000199.2).

Synonyms

HHF5; CD220; Phospho-Insulin Receptor-Y1361

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.