

AP0049

Leader in Biomolecular Solutions for Life Science



Phospho-c-Jun-T239 Rabbit pAb

Catalog No.: AP0049

Basic Information

Observed MW

Refer to Figures

Calculated MW

36kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3725

Swiss Prot

P05412

Immunogen

A synthetic phosphorylated peptide around T239 of human Jun (NP_002219.1).

Synonyms

AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-T239

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.