

AP0305

Leader in Biomolecular Solutions for Life Science



Phospho-AKT2-S474 Rabbit pAb

Catalog No.: AP0305 **1 Publications**

Basic Information

Observed MW

56kDa

Calculated MW

56kDa

Category

Mouse Monoclonal Antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Rat

Background

This gene is a putative oncogene encoding a protein belonging to a subfamily of serine/threonine kinases containing SH2-like (Src homology 2-like) domains, which is involved in signaling pathways. The gene serves as an oncogene in the tumorigenesis of cancer cells. For example, its overexpression contributes to the malignant phenotype of a subset of human ductal pancreatic cancers. The encoded protein is a general protein kinase capable of phosphorylating several known proteins, and has also been implicated in insulin signaling.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

208

Swiss Prot

P31751

Immunogen

A phospho specific peptide corresponding to residues surrounding S474 of human AKT2

Synonyms

PKBB; PRKBB; HIHGHH; PKBBETA; RAC-BETA; Phospho-AKT2-S474

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts from 293 cells using Phospho-AKT2-S474 antibody (AP0305).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% BSA.