

AP0350

Leader in Biomolecular Solutions for Life Science



Phospho-FAM129A-S602 Rabbit pAb

Catalog No.: AP0350

Basic Information

Observed MW

150kDa

Calculated MW

103kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a member of the family with sequence similarity 129 protein family. This gene is highly expressed in several cancer cells and may serve as a prognostic marker for certain cancers. The encoded protein may play a role in regulating p53-mediated apoptosis.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

116496

Swiss Prot

Q9BZQ8

Immunogen

A phospho specific peptide corresponding to residues surrounding S602 of human FAM129A

Synonyms

GIG39; NIBAN; C1orf24; FAM129A; Phospho-FAM129A-S602

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.