

Phospho-IGF1R-Y1135/Y1136 + INSR-Y1150/Y1151 Rabbit pAb

Catalog No.: AP0547

Basic Information

Observed MW

Refer to figures

Calculated MW

154kDa/155kDa/156kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor generates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2014]

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3480/3643

Swiss Prot

P08069/P06213

Immunogen

A synthetic phosphorylated peptide around Y1135 & Y1136 of human IGF1RIGF1R (NP_000866.1).

Synonyms

IGF1R/INSR; Phospho-IGF1R-Y1135/Y1136 + INSR-Y1150/Y1151

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.