

AP0662

Leader in Biomolecular Solutions for Life Science



Phospho-ARAF-S299 Rabbit pAb

Catalog No.: AP0662

Basic Information

Observed MW

Refer to figures

Calculated MW

68kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Enables protein serine/threonine kinase activity. Involved in negative regulation of apoptotic process; regulation of TOR signaling; and regulation of cellular protein metabolic process. Predicted to be active in cytosol and mitochondrion. Biomarker of high grade glioma.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

369

Swiss Prot

P10398

Immunogen

A synthetic phosphorylated peptide around S299 of human ARAF (NP_001645.1).

Synonyms

PKS2; A-RAF; ARAF1; RAFA1; Phospho-ARAF-S299

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.