

AP0855

Leader in Biomolecular Solutions for Life Science



# Phospho-Lyn-Y507 Rabbit pAb

Catalog No.: AP0855

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

59kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human

## Background

---

This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

4067

### Swiss Prot

P07948

### Immunogen

A synthetic phosphorylated peptide around Y507 of human Lyn (NP\_002341.1?).

### Synonyms

JTK8; p53Lyn; p56Lyn; Phospho-Lyn-Y507

## Contact

---

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.