

AP0919

Leader in Biomolecular Solutions for Life Science



Phospho-TAZ-S89 Rabbit pAb

Catalog No.: AP0919

2 Publications

Basic Information

Observed MW

Refer to Figures

Calculated MW

44kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Enables transcription coactivator activity. Involved in several processes, including hippo signaling; positive regulation of cell differentiation; and regulation of signal transduction. Located in cytosol and nuclear body.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

25937

Swiss Prot

Q9GZV5

Immunogen

A synthetic phosphorylated peptide around S89 of human WWTR1 (NP_056287.1).

Synonyms

TAZ; Phospho-TAZ-S89

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.