

AP1153

Leader in Biomolecular Solutions for Life Science



# Phospho-Cdc6-S54 Rabbit mAb

Catalog No.: AP1153 **Recombinant**

## Basic Information

### Observed MW

63kDa

### Calculated MW

63kDa

### Category

SMab Recombinant Monoclonal  
Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Mouse

### CloneNo number

ARC1693

## Background

The protein encoded by this gene is highly similar to *Saccharomyces cerevisiae* Cdc6, a protein essential for the initiation of DNA replication. This protein functions as a regulator at the early steps of DNA replication. It localizes in cell nucleus during cell cycle G1, but translocates to the cytoplasm at the start of S phase. The subcellular translocation of this protein during cell cycle is regulated through its phosphorylation by Cdk. Transcription of this protein was reported to be regulated in response to mitogenic signals through transcriptional control mechanism involving E2F proteins.

## Recommended Dilutions

WB 1:500 - 1:1000

## Immunogen Information

### Gene ID

990

### Swiss Prot

Q99741

### Immunogen

A synthetic phosphorylated peptide around S54 of human Cdc6 (Q99741).

### Synonyms

CDC18L; HsCDC6; MGORS5; HsCDC18; Phospho-Cdc6-S54

## Contact

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

## Validation Data



Western blot analysis of extracts of NIH/3T3 cells, using Phospho-Cdc6-S54 Rabbit mAb (AP1153) at 1:1000 dilution. NIH/3T3 cells were treated by Paclitaxel (1  $\mu$ M) at 37°C for 20 hours. NIH/3T3 cells were treated by Hydroxyurea (4 mM) at 37°C for 20 hours. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 $\mu$ g per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit (RM00020). Exposure time: 5s.