

HRP Goat Anti-Rabbit IgG (H+L)

Catalog No.: AS014 **422 Publications**

Basic Information

Catalog No.

AS014

Observed MW

Calculated MW

Category

Secondary antibody

Applications

ELISA, WB

Cross-Reactivity

Background

Recommended Dilutions

ELISA 1:5000 - 1:10000

WB 1:2000 - 1:10000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Rabbit IgG

Synonyms

Contact

 | www.abclonal.com

Product Information

Source

Goat

Isotype

Horseradish peroxidase
conjugated IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of HeLa cells, using HRP Goat Anti-Rabbit IgG (H+L) antibody (AS014) at 1:4000-1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 3s.