

[One Step]USP25 Antibody Kit

Catalog No.: RK05720

Basic Information

Catalog No.

RK05720

Applications

WB

Cross-Reactivity

Human, Mouse, Rat

Observed MW

130kDa

Calculated MW

122kDa/125kDa/129kDa

Category

Antibody kit

Background

Product Information

Source

Rabbit

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Avoid repeated freeze-thaw cycles.

Contact

 | www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05720-1	USP25 Rabbit pAb	1:1000 dilution
RK05720-2	HRP Goat Anti-Rabbit IgG (H+L)	1:4000 dilution

Immunogen Information

Gene ID

29761

Swiss Prot

Q9UHP3


Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human USP25 (NP_037528.3).

Synonyms

USP25;USP21

Validation Data



Western blot analysis of extracts of various cells, using USP25 antibody (A7975) at 1:1000 dilution ratio through one-step method.