

RK05736

Leader in Biomolecular Solutions for Life Science



[One Step] PGAM1 antibody Kit

Catalog No.: RK05736

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

29kDa

Calculated MW

29kDa

Category

Antibody kit

Background

The protein encoded by this gene is a mutase that catalyzes the reversible reaction of 3-phosphoglycerate (3-PGA) to 2-phosphoglycerate (2-PGA) in the glycolytic pathway. Two transcript variants encoding different isoforms have been found for this gene.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05736-1	PGAM1 Rabbit pAb	29kDa
RK05736-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	29kDa

Immunogen Information

Gene ID

5223

Swiss Prot

P18669

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-254 of human PGAM1 (NP_002620.1).

Synonyms

PGAMA; PGAM-B; HEL-S-35; PGAM1