

RK05737

Leader in Biomolecular Solutions for Life Science



# [One Step] Hsp90 $\beta$ antibody Kit

Catalog No.: RK05737

2 Publications

## Basic Information

### Applications

### Cross-Reactivity

Human, Mouse, Rat

### Observed MW

80kDa

### Calculated MW

83kDa

### Category

Antibody kit

## Background

This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes.

## Product Information

### Source

Rabbit

### Purification

### Storage

Avoid repeated freeze-thaw cycles.

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05737-1	[KO Validated] Hsp90 $\beta$ Rabbit pAb	83kDa
RK05737-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	80kDa

## Immunogen Information

### Gene ID

3326

### Swiss Prot

P08238

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Hsp90 $\beta$  (NP\_031381.2).

### Synonyms

HSP84; HSPC2; HSPCB; D6S182; HSP90B;  $\beta$