

Nuclease-free H₂O



Catalog: RK20102

Size: 25 ml

Concentration: /

Components:

Nuclease-free H ₂ O	RM20214
--------------------------------	---------

Storage: -20°C

QC Process:

- ◆ Free of detectable DNA exonuclease and endonuclease.
- ◆ No residual host genomic DNA is detected by PCR.
- ◆ Free of detectable RNase activities.

Order: order@abclonal.com
Tech: support@abclonal.com

For research purposes only. Not for therapeutic or diagnostic purposes in humans or animals.
Please visit <http://abclonal.com> for a complete listing of recommended products.