

Certificate of Analysis

Product Name:	ABScript III RT Master Mix for qPCR with gDNA Remover				
Cat No.:	RK20429	Lot No.:	962129J18W19	Exp:	2024.04
Conc.:	-			Storage:	-20°C

Assay Name/Specification (Minimum release criteria)	Data	Result
Functional Testing (RT-qPCR) ABScript III RT Mix for qPCR with gDNA Remover is functionally tested in RT-qPCR with human total RNA template, resulting in a standard curve with a calculated qPCR efficiency of 90-110%, and a dynamic range of 5 orders of magnitude.	90%<E<110%	Pass

Conclusion:

In compliance.

Authorized By: Li Shuangjie

Tested By: Wang Yuhong

Date: 2021-10-12

ABclonal Biotechnology Co, Ltd

