

RP00424

Leader in Biomolecular Solutions for Life Science



Recombinant Human Siglec-4a/MAG Protein

Catalog No.: RP00424 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	4099	P20916

Tags

C-His

Synonyms

GMA; S-MAG; SIGLEC-4A; SIGLEC4A; SPG75;MAG;S-MAG;SIGLEC-4A;SIGLEC4A;SPG75

Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of 20 mM PB, 150 mM NaCl, pH7.2. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Contact



www.abclonal.com

Background

This protein is a type I membrane protein and member of the immunoglobulin superfamily. It is thought to be involved in the process of myelination. It is a lectin that binds to sialylated glycoconjugates and mediates certain myelin-neuron cell-cell interactions. Three alternatively spliced transcripts encoding different isoforms have been described for this gene.

Basic Information

Description


Recombinant Human Siglec-4a/MAG Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Gly20-Pro516) of human MAG/Siglec-4a (Accession #P20916) fused with a 6xHis tag at the C-terminus.

Bio-Activity

Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant protein Human MAG/Siglec-4a was determined by SDS-PAGE under reducing conditions with Coomassie Blue, showing a band at 115 kDa.