

# Recombinant Mouse P-Selectin/CD62P Protein

Catalog No.: RP01068 **Recombinant**

## Sequence Information

---

Species	Gene ID	Swiss Prot
Mouse	20344	Q01102

### Tags

C-Fc & 6×His

### Synonyms

P-Selectin;CD62P;SELP;GMP-140

## Product Information

---

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

< 0.1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

## Background

---

### Basic Information

---

#### Description

Recombinant Mouse P-Selectin/CD62P Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Trp42-Ala709) of mouse P-Selectin/CD62P (Accession #NP\_035477.1.) fused with an Fc, 6×His tag at the C-terminus.

#### Bio-Activity

#### Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant Mouse P-Selectin/CD62P  
Protein was determined by SDS-PAGE with  
Coomassie Blue, showing a band at  
130-150kDa.