

A0484

Leader in Biomolecular Solutions for Life Science



Nestin Rabbit pAb

Catalog No.: A0484

5 Publications

Basic Information

Observed MW

280kDa

Calculated MW

177kDa

Category

Polyclonal Antibody

Applications

WB,IF/ICC,FC

Cross-Reactivity

Human

Background

This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells.

Recommended Dilutions

| | |
|---------------|----------------|
| WB | 1:500 - 1:1000 |
| IF/ICC | 1:20 - 1:50 |
| FC | 1:20 - 1:50 |

Immunogen Information

Gene ID

10763

Swiss Prot

P48681

Immunogen

A synthetic peptide of human Nestin

Synonyms

Nbla00170; Nestin

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.

Validation Data



Western blot analysis of various lysates using Nestin Rabbit pAb (A0484) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.