

A0680

Leader in Biomolecular Solutions for Life Science



## PSMG2 Rabbit pAb

Catalog No.: A0680

### Basic Information

**Observed MW**

29kDa

**Calculated MW**

29kDa

**Category**

Mouse Monoclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human, Mouse, Rat

### Background

Enables molecular adaptor activity. Involved in chaperone-mediated protein complex assembly. Located in nucleus. Part of chaperone complex.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

56984

**Swiss Prot**

Q969U7

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human PSMG2 (NP\_064617.2).

**Synonyms**

PAC2; HCCA3; CLAST3; MDS003; PRAAS4; HsT1707; TNFSF5IP1; PSMG2

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

## Validation Data

---



Western blot analysis of extracts of 293T cells, using PSMG2 antibody (A0680).  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.