

A0705

Leader in Biomolecular Solutions for Life Science



CDK6 Rabbit pAb

Catalog No.: A0705

5 Publications

Basic Information

Observed MW

37kDa

Calculated MW

37kDa

Category

Mouse Monoclonal Antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene is a member of the CMGC family of serine/threonine protein kinases. This kinase is a catalytic subunit of the protein kinase complex that is important for cell cycle G1 phase progression and G1/S transition. The activity of this kinase first appears in mid-G1 phase, which is controlled by the regulatory subunits including D-type cyclins and members of INK4 family of CDK inhibitors. This kinase, as well as CDK4, has been shown to phosphorylate, and thus regulate the activity of, tumor suppressor protein Rb. Altered expression of this gene has been observed in multiple human cancers. A mutation in this gene resulting in reduced cell proliferation, and impaired cell motility and polarity, and has been identified in patients with primary microcephaly.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

1021

Swiss Prot

Q00534

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 250-326 of human CDK6 (NP_001250.1).

Synonyms

MCPH12; PLSTIRE; CDK6

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.