

A0844

Leader in Biomolecular Solutions for Life Science



STK11 Rabbit pAb

Catalog No.: A0844

Basic Information

Observed MW

Refer to figures

Calculated MW

49kDa

Category

Mouse Monoclonal Antibody

Applications

WB,IF/ICC,FC

Cross-Reactivity

Human

Background

The protein encoded by this gene is a serine/threonine kinase that regulates cell polarity and energy metabolism and functions as a tumor suppressor. Mutations in this gene have been associated with the autosomal dominant Peutz-Jeghers syndrome, as well as with skin, pancreatic, and testicular cancers.

Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:100
FC	1:20 - 1:50

Immunogen Information

Gene ID

6794

Swiss Prot

Q15831

Immunogen

A synthetic Peptide of human STK11

Synonyms

PJS; LKB1; hLKB1; STK11

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.