

A10003

Leader in Biomolecular Solutions for Life Science



# EV71 3C Rabbit pAb

Catalog No.: A10003 **4 Publications**

## Basic Information

---

### Observed MW

Refer to Figures

### Calculated MW

### Category

Polyclonal Antibody

### Applications

WB, ELISA

### Cross-Reactivity

Enterovirus

## Background

---

## Recommended Dilutions

---

WB 1:1000 - 1:5000

## Immunogen Information

---

Gene ID Swiss Prot

### Immunogen

Recombinant protein of Enterovirus EV71 3C

### Synonyms

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.