

A10004

Leader in Biomolecular Solutions for Life Science



ALKBH6 Rabbit pAb

Catalog No.: A10004

Basic Information

Observed MW

26kDa

Calculated MW

26kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable dioxygenase activity and metal ion binding activity. Located in focal adhesion and nucleoplasm.

Recommended Dilutions

WB 1:200 - 1:1000

Immunogen Information

Gene ID

84964

Swiss Prot

Q3KRA9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 97-266 of human ALKBH6 (NP_116267.3).

Synonyms

ABH6; ALKBH6

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using ALKBH6 antibody (A10004) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.