

LY75 Rabbit pAb

Catalog No.: A10065

Basic Information

Observed MW

250kDa

Calculated MW

198kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse

Background

Predicted to enable signaling receptor activity. Predicted to be involved in immune response and inflammatory response. Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

4065

Swiss Prot

O60449

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 216-501 of human LY75 (NP_002340.2).

Synonyms

CD205; LY-75; CLEC13B; DEC-205; GP200-MR6; LY75

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates, using LY75 antibody (A10065) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 180s.