

A10083

Leader in Biomolecular Solutions for Life Science



INSRR Rabbit pAb

Catalog No.: A10083

Basic Information

Observed MW

144kDa

Calculated MW

144kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Enables transmembrane receptor protein tyrosine kinase activity. Involved in actin cytoskeleton reorganization; cellular response to alkaline pH; and protein autophosphorylation. Part of receptor complex.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

3645

Swiss Prot

P14616

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 747-921 of human INSRR (NP_055030.1).

Synonyms

IRR; INSRR

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using INSRR Rabbit pAb (A10083) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates / proteins: 25 µg per lane.
Blocking buffer: 3 % nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.