

A1015

Leader in Biomolecular Solutions for Life Science



Alpha-1 Antitrypsin (SERPINA1) Rabbit pAb

Catalog No.: A1015

4 Publications

Basic Information

Observed MW

51kDa

Calculated MW

47kDa

Category

Polyclonal Antibody

Applications

WB, IHC-P, IF/ICC, ELISA

Cross-Reactivity

Human, Mouse

Background

The protein encoded by this gene is a serine protease inhibitor belonging to the serpin superfamily whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. This protein is produced in the liver, the bone marrow, by lymphocytic and monocytic cells in lymphoid tissue, and by the Paneth cells of the gut. Defects in this gene are associated with chronic obstructive pulmonary disease, emphysema, and chronic liver disease. Several transcript variants encoding the same protein have been found for this gene.

Recommended Dilutions

| | |
|--------|---------------|
| WB | 1:100 - 1:500 |
| IHC-P | 1:50 - 1:200 |
| IF/ICC | 1:50 - 1:200 |

Immunogen Information

Gene ID

5265

Swiss Prot

P01009

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 156-255 of human Alpha-1 Antitrypsin (SERPINA1) (NP_000286.3).

Synonyms

PI; A1A; AAT; PI1; A1AT; nNIF; PRO2275; alpha1AT; Alpha-1 Antitrypsin (SERPINA1)

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.


Validation Data



Western blot analysis of extracts of various cell lines, using Alpha-1 Antitrypsin (SERPINA1) antibody (A1015) at 1:500 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 60s.



Immunohistochemistry analysis of paraffin-embedded mouse kidney using Alpha-1 Antitrypsin (SERPINA1) Rabbit pAb (A1015) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunofluorescence analysis of HepG2 cells using Alpha-1 Antitrypsin (SERPINA1) Rabbit pAb (A1015) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.