

A10173

Leader in Biomolecular Solutions for Life Science



## PDX1 Rabbit pAb

Catalog No.: A10173

4 Publications

### Basic Information

#### Observed MW

39kDa

#### Calculated MW

31kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human, Mouse

### Background

The protein encoded by this gene is a transcriptional activator of several genes, including insulin, somatostatin, glucokinase, islet amyloid polypeptide, and glucose transporter type 2. The encoded nuclear protein is involved in the early development of the pancreas and plays a major role in glucose-dependent regulation of insulin gene expression. Defects in this gene are a cause of pancreatic agenesis, which can lead to early-onset insulin-dependent diabetes mellitus (IDDM), as well as maturity onset diabetes of the young type 4 (MODY4).

### Recommended Dilutions

WB 1:500 - 1:1000

### Immunogen Information

#### Gene ID

3651

#### Swiss Prot

P52945

#### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PDX1 (NP\_000200.1).

#### Synonyms

GSF; IPF1; IUF1; IDX-1; MODY4; PDX-1; STF-1; PAGEN1; PDX1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification


Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

# Validation Data



Western blot analysis of extracts of various cell lines, using PDX1 antibody (A10173) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 90s.