

A10301

Leader in Biomolecular Solutions for Life Science



NRBP1 Rabbit pAb

Catalog No.: A10301

Basic Information

Observed MW

60kDa

Calculated MW

60kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to enable protein homodimerization activity. Involved in endoplasmic reticulum to Golgi vesicle-mediated transport. Located in endomembrane system.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

29959

Swiss Prot

Q9UHY1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 406-535 of human NRBP1 (NP_037524.1).

Synonyms

MADM; NRBP; BCON3; MUDPNP; NRBP1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using NRBP1 antibody (A10301) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.