

A10424

Leader in Biomolecular Solutions for Life Science



ENTPD3 Rabbit pAb

Catalog No.: A10424

Basic Information

Observed MW

59kDa

Calculated MW

59kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

This gene encodes a plasma membrane-bound divalent cation-dependent E-type nucleotidase. The encoded protein is involved in the regulation of extracellular levels of ATP by hydrolysis of it and other nucleotides. Multiple transcript variants have been described.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

956

Swiss Prot

O75355

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 45-366 of human ENTPD3 (NP_001278890.1).

Synonyms

HB6; CD39L3; NTPDase-3; ENTPD3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of lysates from mouse lung, using ENTPD3 Rabbit pAb (A10424) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.