

# RELL2 Rabbit pAb

Catalog No.: A10620

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

32kDa

### Category

Primary antibody

### Applications

ELISA

### Cross-Reactivity

Human

## Background

Predicted to enable collagen binding activity. Involved in positive regulation of p38MAPK cascade. Predicted to be located in extracellular matrix.

## Recommended Dilutions

## Immunogen Information

### Gene ID

285613

### Swiss Prot

Q8NC24

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 44-303 of human RELL2 (NP\_776189.3).

### Synonyms

C5orf16; RELL2

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.