

A1072

Leader in Biomolecular Solutions for Life Science



## RNASE12 Rabbit pAb

Catalog No.: A1072

### Basic Information

#### Observed MW

#### Calculated MW

17kDa

#### Category

Polyclonal Antibody

#### Applications

ELISA

#### Cross-Reactivity

Human

### Background

Predicted to enable nucleic acid binding activity. Predicted to be located in extracellular region.

### Recommended Dilutions

### Immunogen Information

#### Gene ID

493901

#### Swiss Prot

Q5GAN4

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human RNASE12 (NP\_001019993.1).

#### Synonyms

RAI1; HEL-S-85p; RNASE12

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.